

T-TAPE COMPANY

INTRODUCTION



T-TAPE COMPANY

Introduction:

In the realm of material science and manufacturing, innovation is the driving force behind progress and growth. The T Tape Company (located in Putte, the Netherlands), a prominent manufacturer of the low temperature thermoplastic materials, stands as an example of how technological advancements can revolutionize the medical industry. Through its unwavering commitment to research, development, and production, T Tape Company has paved the way for the adoption of low temperature thermoplastics in various healthcare sectors, redefining the possibilities of material applications and transforming medical care worldwide



T Tape Company BV:

- ▶ T Tape is a privately-owned company established in 1987.
- ▶ It manufactures high-end products developed in close cooperation with medical professionals in the fields of hand therapy, rehabilitation, orthopedics, and radiotherapy.
- ▶ The management team has many years of technical and clinical background.
- ▶ The founder has a proven track record in the development of innovative thermoplastic materials for more than 40 years.
- ▶ The company follows a lean business model, focusing on product development, marketing, sales, quality assurance & regulatory compliance, and logistics.
- ▶ Sales are conducted through international distributors, both direct (in 35+ countries) and indirect (in 60+ countries).



T-TAPE COMPANY



T-TAPE COMPANY


The Genesis of LTTPS

- ▶ Low temperature thermoplastic materials (LTTTP) are a class of polymers that exhibit unique properties when heated to relatively low temperatures, allowing them to be molded, shaped, and reshaped multiple times without undergoing degradation. The origins of these materials can be traced back to the mid-20th century when researchers began to explore the potential of polymers and their behavior under varying temperature conditions.
- ▶ T Tape Company emerged as a pioneering force in this field, leveraging extensive research and development efforts to perfect the formulation and production of low temperature thermoplastics. By refining the polymer composition and optimizing manufacturing processes, the company successfully harnessed the potential of these materials, setting the stage for a revolutionary leap in material science.



MODERN LOW-TEMPERATURE THERMOPLASTIC MATERIALS

- ▶ **Turbocast®**
- ▶ **Tubocast Ortho (NS)**
- ▶ **IMMO+ (NS)**
- ▶ **Beachcast**
- ▶ **B-Cast**



T-TAPE COMPANY

TWO BUSINESS UNITS


**RADIOTHERAPY &
NUCLEAR MEDICINE**



ORTHOPAEDICS





T-TAPE COMPANY



THERMOPLASTICS FOR IMMOBILIZATION IN RADIOTHERAPY:

Required to establish and maintain the patient in a fixed, well-defined position over the course of multiple radiation treatments, or to prevent the patient from moving during a single treatment session




T-TAPE COMPANY

THERMOPLASTICS FOR IMMOBILIZATION IN ORTHOPAEDICS:

fracture immobilization
static and dynamic splinting
splinting for burned patients
post-operative splinting
protective devices




LOW TEMPERATURE **THERMOPLASTICS**



GENERAL CHARACTERISTICS

- ▶ Does not emit toxic or noxious gases or vapors during application
- ▶ X-ray transparent
- ▶ After application, the material transforms into a hard splint that does not require reinforcement
- ▶ Due to its exceptional properties, it follows the exact contours of the treated surface
- ▶ The finished device is lightweight
- ▶ Programmed plastic memory allows the product to return to its original shape for further adjustments



TURBOCAST®

Advantages and characteristics

ADVANTAGES OF THE **FOAM** COATING

- ▶ Turbocast® does not adhere to itself in hot water, allowing sizable surfaces to be placed in a small hot-water container.
- ▶ It only sticks when firmly pressed together, preventing accidental bonding.
- ▶ The foam coating creates a thin layer between the patient's skin and the thermoplastic, reducing the risk of skin maceration.



MOLDING CHARACTERISTICS



| Transparency | Memory | Drape | Stretch | Ease of finishing | Bounding |
|--------------|--------|---------|----------|-------------------|---|
| Moderate | 100% | Maximum | Moderate | Maximum | Moderate, no possibility of accidental bounding |

TURBOCAST® IS
MANUFACTURED
IN THE FOLLOWING
COLORS




- ▶ Peach (Skin)
- ▶ Yellow
- ▶ Green
- ▶ Blue
- ▶ Fuchsia-red




TURBOCAST® IS
AVAILABLE
IN THE FOLLOWING
PERFORATIONS




MICRO
0,8 -1,6 - 2,0 mm




OPTI
2,0 - 2,5 - 3,0 mm



MINI
2,0 - 2,5 - 3,0 - 4,0
mm



MULTI
2,0 - 2,5 - 3,0 mm



TURBOCAST® ORTHO (NS)

Advantages and characteristics

TURBOCAST® ORTHO


- ▶ Transparent-when-heated material.
- ▶ Easy auto-bounding (no solvents needed).
- ▶ Lower processing temperature (65°C - 159°F).






MOLDING CHARACTERISTICS

| Transparency | Memory | Drape | Stretch | Ease of finishing | Bounding |
|--------------|--------|----------|---------|-------------------|----------|
| Maximum | 100% | Moderate | Maximum | Maximum | Maximum |



IMMO+ (NS)

Advantages and characteristics



Special advantages

- Extra rigidity
- Extra strength
- Non sticky surface



Colours

IMMO+ (NS) is manufactured in WHITE, RED,
ROYAL BLUE and BLACK





MOLDING CHARACTERISTICS

| Transparency | Memory | Drape | Stretch | Bounding |
|--------------|--------|---------|---------|----------|
| Maximum | 100% | Maximum | Maximum | Minimum |



PROCESSING METHODS

- ▶ Water (from 65°C)
- ▶ Heat gun
- ▶ Heating plate / table

PROCESSING METHODS

The material can be processed after 2-minutes heating at 65°C or 1-minute heating at 75°C or until a dark “glazy” discoloration has appeared.



INSTRUCTIONS OF USE

After heating in water at 65°C lift the material out and place it immediately on a dry towel. Fold half of the towel over the whole plate to dab excess moisture.



INSTRUCTIONS OF **USE**

The material can now be placed directly on bare skin. By this time, the temperature will have dropped to approximately 40°C/105°F, which is generally considered comfortable.



INSTRUCTIONS OF **USE**

If working at a higher temperature is required in order to get more modelling time, it is advisable to protect the skin with stockinet or another kind of tubular bandage.



USING A **WATER BATH** (TURBOCAST® ORTHO)

- ▶ To soften Turbocast® Ortho, immerse the material in a water bath at a temperature of 60°C - 65°C (140°F - 149°F) for 60-90 seconds or until it reaches the desired softness.
- ▶ Once ready, remove the material from the water and place it on a flat surface, preferably a wooden board, for approximately 10 seconds. The temperature will quickly decrease to approximately 35°C (95°F).



T-TAPE COMPANY



THANK YOU!

You can contact me at turbocast@planet.nl